HOW to STORE and MAINTAIN NOVASINA HUMIDITY STANDARDS

Novasina is supplier of saturated pure chemical salt solutions, which can be used over many years as often as desired to check and calibrate Novasina instruments, as long as these salts are carefully maintained.

A humidity standard is correctly saturated, when there is a visible surplus of water beside the salt crystals (only exception: salt 11% may look dry), and the salt should be „slushy“. Then the relative humidity value is correct.

Storage:

We recommend to store the salts in a cabinet as used for light chemicals. The salts of different humidity values should not be kept together in a common container, otherwise they try to equilibrate between themselves. The salt values below the ambient air humidity, mainly 6%, 11%, 33%, will absorb water. Salts above the ambient air humidity, as 90% and 98%, lose water and will dry out.

All Novasina salts are supplied in a well sealing lab-quality container, thus avoiding a humidity exchange with the ambient air for a very long period of time.

After use, the salts should immediately be placed back into the container.

Tropical climate:

we recommend to store low humidity salts in their container in a dessicator, perhaps together with silicagel.

Humidity check on a HygroMate handheld instrument

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Novasina SC Humidity check for industrial air application.

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
Novasina SAL-T humidity check for water activity application.

Saturate dry salts again as follows:

- Sensor Checks SC in tubular form can be opened (screw on top) and distilled water injected with a syringe, until a surplus of water is visible.
- salt tablets SAL-T: use a syringe to inject distilled water through the small hole at the rim (sealed with white silicone paste). Then seal the hole again well (we recommend Dow Corning 734).

Both humidity standards are comfortable to use, because there is no direct contact with the salt. The water molecules wander through the special membrane, but droplets are stopped.

Low value salts, mainly 6% and 33%, absorb more and more water which eventually dilutes the salt, until all crystals are *dissolved*. In this case, the humidity standard is *useless*. It must be replaced by a new one, before all salt crystals disappear!

For accurate results, the salt shall not form a solid block. Therefore **ALWAYS SHAKE** each salt very well from time to time, especially the 98%, but latest before calibrating an instrument.

Please also refer to the safety instructions delivered with each salt.

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