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
novasina
The Art of Precision Measurement

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Novasina humidity standards have several critical advantages over unsaturated salt solutions.

Novasina humidity standards are totally hassle-free and can be used time and time again. For maximum service life, the humidity standards should be hermetically resealed immediately after use. If they do eventually dry out, all they need is to be moistened again with distilled water until visibly saturated. A drop or two of water in the humidity standards is an indication that the humidity reading is accurate.

Compare these with unsaturated humidity standards like lithium chloride. They are intended for once-only use and it is impossible to tell whether they are still useable or to see the nominal humidity value simply by looking at them. It also means that every calibration creates toxic waste that requires special disposal.

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1 & 2 Giftig/Toxic/Toxique
3 & 4 Gesundheitsschädlich/Injurious to health/Nocif

²⁾ Literatur-Referenz /Literature reference/Référence littérature

- A: Greenspan, Humidity Fixed Points of Binary Saturated Aqueous Solutions. Journal of Research of the National Bureau of Standards. Vol. 81A, No.1, Jan.–Febr. 1977
B: Robinson R.A. and Stokes R.H. Electrolyte Solutions, Butterworths London 1959
C: Wexler, Humidity and Moisture, Vol. 3, Fundamentals and Standards, Robert E. Krieger, Publishing Company, Huntington, New York

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Safe disposal of Novasina humidity standards:

Unused product: Dispose of as special waste as stipulated by local and national legal regulations.

Waste code 2640 (VVS-CH)
Abfallschlüssel no. 51523 (Germany).

Non-detoxified packaging:
Disposal as for unused product.

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Manufacture : Novasina AG