



**Accessories for LabSwift-aw, ms-1 aw
and LabMaster/PARTNER-aw**



111 9979 SI Set Master

SI, sorption isotherm set


for the LabMASTER-aw and LabPARTNER-aw

This accessory, together with the Novasina humidity standards SAL-T, makes it possible to easily determine several points of a sorption isotherm. Incl. a plastic case for transportation.

In addition, a precision balance and an oven are required.

**SI-Set MASTER/PARTNER
complete with:**

- 1 large cup insert
- 1 substance basket
- 1 weighing glass
- 1 pair of tongs



- 111 1001** -> eVC-21
- 111 1003** -> eVC-26
- 111 7212** -> Redox
- 111 0995** -> eVALC
- 111 0997** -> eVMT-2
- 111 0999** -> eVMS-2

Cell protection filters

for all LabMASTER-aw and LabPARTNER-aw.

A number of filters is available to protect the measurement cell from aggressive vapours or abrasion.

-> Please consult the appropriate datasheets

Filters for enBSK head:

- 111 1001** -> eVC-21 (10 pieces)
- 111 1003** -> eVC-26 (10 pieces)
- 111 7212** -> Redox (1 piece)
- 111 0995** -> eVALC (2 pieces)
- 111 0997** -> eVMT-2 (5 pieces)
- 111 0999** -> eVMS-2 (5 pieces)

Weight: 30g



111 0601 ePW sample dishes

Sample dishes ePW

For all LabMASTER-aw and LabPARTNER-aw

Box with 100 sample cups.

These cups are for all Novasina aw instruments. They are disposable, intended for a single use.

A box contains:

100 sample cups ePW

Material : Polypropylene

Dimensions : 40mm diameter
x 12mm deep