

World New Innovation

Water activity measurement instrument

NOVASINA is the leader since more than 30 YEARS of precision **WATER ACTIVITY** measurement systems in foodstuff as well

as drugs, chemical products and cosmetics. Best "Swiss quality" thanks experience, research and development

a_w



Advantages:

- ↪ Precise measurement
- ↪ Worldwide unique "Novalyte" measurement technology
- ↪ Highest reliability and long term stability
- ↪ Simple to use and handle
- ↪ Intuitive menu structure
- ↪ High precision temperature controll chamber
- ↪ Built in "pre-conditioning chamber"
- ↪ Large-scale display, backlighted
- ↪ 7 point Factory calibration
- ↪ Checking, testing and re-adjusting with SAL-T humidity standards
- ↪ Modular design
- ↪ Software updates possible at any time
- ↪ Easy on site service possibilities

Novasina AG is proud to present its newest product for water activity measurement.

The **LabMASTER-aw** and **LabPARTNER-aw** product series impress with their unique accuracy, reproducibility, robust-ness and modularity.

With that the new system is setting new standards and fulfils already today the regulations and recommendations of international commissions regarding trace ability and safety.

The **LabMASTER-aw** is available in three versions and can control up to **9 LabPARTNER-aw** by a electronic bus-system. Thanks to a clever software design updates can be easily performed at any time.



novasina
The Art of Precision Measurement

Neuheimstrasse 12
CH-8853 Lachen SZ
Phone: +41 (0)55 642 67 67
Fax: +41 (0)55 642 67 70
Mail: info@novasina.ch

www.novasina.ch